


QUICK START GUIDE

CSC Scientific Economy Moisture Balance




How to Start Drying

<ol style="list-style-type: none"> 1. Put an empty pan on the balance Close the lid 2. Press Start Stop 3. Press →0/T+ to zero the balance 	<ol style="list-style-type: none"> 5. Open the lid 6. Place the sample on the pan 7. Close then lid and the test will automatically begin
---	--

How to Set the Temperature

<ol style="list-style-type: none"> 1. Press 2. Press one time 3. Press two times 	<ol style="list-style-type: none"> 5. When finished press three times 6. Then press or to select the desired temperature
--	---

How to Set the Auto Finish




(Amount of time before the test shuts off after the weight stops changing 1mg or less)

SETTINGS: 1= 20sec 2= 50sec 3= 120sec 4= 180sec 5= 240sec Manual= Timed=

<ol style="list-style-type: none"> 1. Press 2. Press one time 	<ol style="list-style-type: none"> 3. Then press or to select the desired setting 4. When finished press two times
---	---

How to Change the Read Out Formula



$\frac{0}{0}^M$ = % Moisture [(Start Weight-End Weight)/Start Weight]

$\frac{0}{0}^D$ = % Solids [End Weight/Start Weight]

$\frac{0}{0}^R$ = % Dry Weight Basis [(Start Weight-End Weight)/End Weight]

<ol style="list-style-type: none"> 1. Press Test menu 2. Press two times 3. Press one time 	<ol style="list-style-type: none"> 4. Then press or to select the desired setting 5. When finished press two times
--	---

How to Switch to a Timed Test

<ol style="list-style-type: none"> 1. Press Test menu 2. Press one time 3. Press or until a picture of a clock is displayed on the LED 4. Press one time 	<ol style="list-style-type: none"> 5. Press or to adjust the hours 6. Press one time to adjust the minutes 7. When finished press three times
---	--

How to Set Up the Printer

1. Plug the Printer into the “RS 232 Port” on the back of the moisture balance.

2. The unit prints automatically but you may also press at any time to print.

CSC Scientific Company, Inc. Phone: 1-800-621-4778 Fax: 1-703-280-5142

Email: info@cscscientific.com web: <http://www.cscscientific.com>